EN

TotalCare System

TotalCare System is a simple and convenient way to care for your hard and rigid gas permeable contact lenses.
The **TotalCare** System is formulated to help increase comfort and helps maximize lens wearing time.
The **TotalCare** System consists of the following components:

- TotalCare** Daily Cleaner
- TotalCare** Disinfecting/Storing/Wetting Solution
- TotalCare** Lens Case (non-sterile)

Daily Care Routine

Always wash and dry your hands thoroughly before handling your contact lenses.

Cleaning

- Remove the left lens from your eye and place on the palm of your hand.
- Place 2-3 drops of **TotalCare** Daily Cleaner on each lens surface and rub for 15-20 seconds with your finger. Never rub in a circular motion as it may not effectively clean the lens edge.
- Rinse each side of the lens using **TotalCare** Disinfecting, Storing and Wetting Solution or sterile saline solution.
- Place the lens in the holder of the left cap (L) and close the lens case with this cap.
- Fill the lens case with **TotalCare** Disinfecting, Storing and Wetting Solution. Ensure the lens case is completely filled to the top.
- Repeat this procedure with the right contact lens.

Disinfecting

- Leave your lenses to soak in **TotalCare** Disinfecting, Storing and Wetting Solution for at least 6 hours or overnight to ensure thorough disinfection. Store the lenses in the unopened lens case until ready to wear.
- Rinse your lenses with **TotalCare** Disinfecting, Storing and Wetting Solution.
- Your lenses are now ready to wear.
- Discard the used solution, rinse your lens case with **TotalCare** Disinfecting, Storing and Wetting Solution and leave to air dry.

Storage

- Lenses can be stored for up to 30 days in the unopened lens case.
- If they are stored for longer periods of time, disinfect again before inserting your lenses.

Wetting

- To wet your contact lenses place 1-2 drops of **TotalCare** Disinfecting, Storing and Wetting Solution on the inner surface of each contact lens before inserting.

Protein Removal

- The proteins that occur naturally in your tear film can build up a deposit on your contact lenses during normal wear. These deposits can reduce wearing comfort, vision and the lifespan of your lenses.
- Where required, the weekly use of protein remover tablets provides optimal care for your lenses, please contact your eye care professional for advice.

Precautions

All instructions must be followed. Failure to do so may result in serious eye infection or injury.

General

Do not use any of these products if you are allergic to any of their contents.
If eye irritation occurs discontinue lens wear and consult your eye care practitioner.
Do not use if tamper-evident seal is broken or missing.

Use

The **TotalCare** System is not intended for use with soft contact lenses.

Always use fresh solution.
Do not use after the expiry date.

Do not insert lens directly into the eye after soaking the lens in solution containing a protein remover tablet.
To avoid contamination, do not touch bottle tip on any surface and avoid contact with eyes.
Do not swallow the tablets or solutions.

Do not rinse your lenses or lens case with tap water.

It is recommended you use **TotalCare** Disinfecting, Storing and Wetting Solution with the **TotalCare** lens case.

Do not change contact lens care system without consulting an eye care professional.

Storage

Store all products in the **TotalCare** System in a dry place.

Store between 5°C and 25°C.

Close the bottles tightly immediately after use.

Keep out of reach of children.

Discard

Discard any remaining **TotalCare** Disinfecting, Storing and Wetting Solution 45 days after first opening the bottle.

Discard any remaining **TotalCare** Daily Cleaner 60 days after first opening the bottle.

Contraindications

None known.

Product:

TotalCare Disinfecting,
Storing and Wetting Solution:

Contains:

Polyhexamethylene biguanide (PHMB) - 0.0005%
Edetate disodium
Hydroxyethyl cellulose
Sodium phosphate
Sodium chloride
Hydrochloric acid
Purified water

Polyhexamethylene biguanide - 0.0001%
Tyloxapol - 0.025%
Hydroxypropyl methylcellulose - 0.15%
Tromethamine
Edetate disodium
Hydrochloric acid
Sodium chloride
Purified water

Note: It is possible that not all products listed above are contained in this pack but are available separately, if required. Please refer to the instructions for use contained in each of these products.

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© Johnson & Johnson Surgical Vision, Inc. 2024