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AVIZOR CLEANING, DISINFECTING, RINSING, WETTING, AND
gp system STORING FOR ALL TYPE OF RIGID AND GAS
 PERMEABLE CONTACT LENSES.

INSTRUCTIONS FOR USE: 1. Wash your hands carefully with water and soap before handling your lenses. 2. Place the lens on the palm of your hand. 3. Add a few drops of AVIZOR GP CLEANER SOLUTION on the surface of each lens. 4. Rub gently with the tip of your finger until a little foam appears. 5. Rinse your lenses abundantly with AVIZOR GP CONDITIONER SOLUTION or AVIZOR SALINE SOLUTION to eliminate all residue of the cleaning solution. 6. Place your lenses in the lens case and cover them with AVIZOR GP CONDITIONER SOLUTION. 7. To guarantee an adequate disinfection, lenses should be soaked in that solution for at least 4 hours or overnight. 8. Before fitting your lenses, add a few drops of AVIZOR GP CONDITIONER on each lens. 9. Keep your lens case clean, and use a fresh supply of solution every day. Last revision date: 2009/11.

PRECAUTIONS: • AVIZOR GP CLEANER SOLUTION causes irritation if not rinsed abundantly. The solution should never come in contact with your eye. • Once opened, the bottles should not be used after 60 days. • Keep bottles well closed after use and keep away from children's reach. • You may store your lenses up to 30 days in the closed lenscase. • Do not touch the dropper or dirty it. • In case of eye irritation consult a specialist. • Do not use these products with soft contact lenses. • Do not use the product if the safety seal is missing or broken. • Never use tap water to rinse your lenses or lens case. • Do not use after the expiration date. • Keep at room temperature. • Before using any ophthalmic medicine, inform your doctor you are a contact lens user. • Do not use if you are allergic to any of the components. • Discard the solution in the lenscase after each use.

COMPOSITION: AVIZOR GP CLEANER: Anionic detergent 10% in aqueous solution. AVIZOR GP CONDITIONER: Viscosity agent, EDTA 0,1%, Polyhexanide 0,0002% in an isotonic, buffered and sterile solution.

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