AVIZOR

lipid
clean

Lipids Cleaning

Specialty for hydrogel
silicone contact lenses

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gb

DAYLY CLEANING SOLUTION FOR
ALL TYPES OF SOFT CONTACT
LENSES

INSTRUCTIONS FOR USE:

1. Wash your hands carefully with water and soap before handling your lenses.
2. Place the lens on the palm of your hand.
3. Shake the bottles and add a few drops of **Avizor Lipid Clean** on the surface of each lens.
4. Rub gently with the tip of your finger.
5. Rinse your lenses abundantly with **Avizor Salina** solution to eliminate all residue of the cleaning solution.

Last revision date: 2009/11

PRECAUTIONS:

- The solution should never come in contact with your eye.
- Once opened, the bottle should not be used after 60 days.
- Maintain the bottle well closed after use.
- Do not use the product if the safety seal is missing or broken.
- Do not ingest the solution and keep away from children's reach.
- Do not allow the dropper to come into contact with any surface to avoid contamination.
- Do not use after the expiration date.
- Keep at room temperature.
- Before using any ophthalmic medicine, inform your doctor you are a contact lens user.

COMPOSITION:

Orgasol 2%, Poloxamer 0.5%, EDTA 0.1% and Polyhexanide 0.0005% in a buffered and isotonic solution.

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AVIZOR S.A.